

CAPP

D E N M A R K

Capp *Maestro*

ELECTRONIC PIPETTES

Capp **Maestro** is the electronic pipette that combines all benefits that the modern researcher needs. Users of electronic pipettes want to have the ability to work as fast as on a mechanical pipette, yet benefiting from the functions of an electronic pipette.

The Favorites button is the first function that you will see on Capp **Maestro's** fully colored OLED display. The colored OLED display gives better optical control while pipetting as they have better contrast and greater screen refresh rate compared to the other LCD screens. Also, OLED display consume less power compared to the other display.

Lightweight construction and elegant ergonomic design with wide hook allows you to rest the pipette at the palm of your hand with no need to actually "grab it".

Capp **Maestro's** tip ejector also works as a finger rest for your thumb. Simply rest you thumb on the tip ejector, and gently activate the pipetting button at the same time.



PIPETTING



CAPP APS
 LANDBRUGSVEJ 10
 DK-5260 ODENSE
 DENMARK

TEL. (+45) 6613 6140
 FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK
 WEB: WWW.CAPP.DK

Capp *Maestro*

ELECTRONIC PIPETTES



Optimize your serial work with advanced functions such as multi dispensing, dilution and mixing, all available on Capp *Maestro* while easily navigating through the user-friendly graphic intuitive interface that is self explanatory and can be operated very easily!

For your routine lab work, Capp *Maestro* allows you to easily store your favorite protocols and allows you to switch from one to the other in a matter of 2 clicks.

If you are tired of exchanging the batteries too often, then try Capp *Maestro* that uses lithium batteries which are ideal for long lasting performance. You can also use this electronic pipette while on charging.

Calibration and password setting functions are also available on this pipette.

Ordering Info

Cat No.	Description	Innaccuracy %	Imprecision %
M10-1	Capp Maestro single channel, 0,2-10ul	2.50 / 1.00	1.80 / 0.50
M100-1	Capp Maestro single channel, 5-100ul	2.50 / 0.70	1.00 / 0.20
M200-1	Capp Maestro single channel, 10-200ul	2.50 / 0.60	0.70 / 0.20
M1000-1	Capp Maestro single channel, 50-1000ul	2.00 / 0.80	0.60 / 0.20
M5000-1	Capp Maestro single channel, 500-5000ul	1.20 / 0.60	0.30 / 0.15

PIPETTING