

# Capp Rondo

CLINICAL CENTRIFUGE

**Capp Rondo Clinical Centrifuge** is a compact and reliable low-speed centrifuge, ideal for blood centrifuge applications involving volumes below 15mL. It stands out in its category for having great features that allow a reliable performance, ease of use and safety.

#### Ease of use:

Rotor capacity of up to 6 x 15mL tubes on a fixed angle rotor.  
Adaptors included to allow usage with most standard collection tubes: 5mL, 7mL, 10ml and 1.5mL/2.0mL.  
Large digital display. Timer setting of 0.5 to 30 minutes.  
Only one click to change from RPM to RCF.

#### Reliable performance:

Smart air flow design reduces noise.  
Precise RPM setting with increments from 100 RPM and setting from 500-6500 RPM.  
Built with a long-life DC motor and a microprocessor.

#### Safety:

Automatic imbalance detection.  
Lid lock safety: lid opens automatically after the operation is completed.



#### Specifications:

Maximum speed RPM / RCF (g)	6500 RPM / 3684 g
Speed setting	500 to 6500 rpm (increments of 100 rpm)
Timer setting	0,5 to 30 minutes or infinite mode
Imbalance detection	Yes
Dimensions	260 x 244 x 205 mm
Noise level	< 60 dB
Weight	4 kg
Rotors and adaptors included	6x 10/15mL, adaptors for 5mL, 7mL, 1.5mL and 2mL tubes

## Ordering Info

Cat No.	Description
CRC-656	Capp Rondo Clinical Centrifuge 6500rpm / 3684g for 6x 15mL