

# CAPP

D E N M A R K

## Capp *Bravo*

SINGLE CHANNEL PIPETTES

Manual pipetting has never been easier! Capp Denmark has combined lightweight ergonomic design with the latest innovations in manual liquid handling to create Capp *Bravo*, one of the most user friendly manual pipettes on the market today.

Capp *Bravo* is constructed of the highest quality autoclavable materials and features feather light plunger and ejection forces. There is a built-in oversized indicator for easy volume monitoring and precise adjustment. With Capp *Bravo* all of these features are packaged in a modern, streamlined, world-renowned Danish design. Experience Capp's symphony of precise measurements and see for yourself why it is called *Bravo!*


**Safety:**  
Fully Autoclavable  
Positive indicator lock

**Ergonomics:**  
Ergonomically constructed hook for comfortable handling  
Smooth tip mounting and ejection

**Controls:**  
Oversized indicator  
Reduced force feather touch plunger and ejector

**Unique Features:**  
Metal shaft with double o-rings for ultra micro sizes  
Ultra lightweight construction



 CAPP APS  
 LANDBRUGSVEJ 10  
 DK-5260 ODENSE  
 DENMARK

TEL. (+45) 6613 6140  
 FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK  
 WEB: WWW.CAPP.DK



## Ordering Info

### Capp *Bravo* variable volume

Cat. No.	Description	Inaccuracy %	Imprecision %	Color
B02-1	Capp pipette, variable vol. 0.1-2 ul, Bravo	12.00 / 2.00	6.00 / 1.00	
B10-1	Capp pipette, variable vol. 0.5-10 ul, Bravo	2.50 / 1.00	1.80 / 0.50	
B20-1U	Capp pipette, variable vol. 2-20 ul, for use with 10ul tips, Bravo	5.00 / 1.00	1.50 / 0.30	
B20-1	Capp pipette, variable vol. 2-20 ul, Bravo	5.00 / 1.00	1.50 / 0.30	
B50-1	Capp pipette, variable vol. 5-50 ul, Bravo	2.00 / 1.00	1.20 / 0.30	
B100-1	Capp pipette, variable vol. 10-100 ul, Bravo	2.50 / 0.70	1.00 / 0.20	
B200-1	Capp pipette, variable vol. 20-200 ul, Bravo	2.50 / 0.60	0.70 / 0.20	
B1000-1	Capp pipette, variable vol. 100-1000 ul, Bravo	2.00 / 0.80	0.60 / 0.20	
B5000-1	Capp pipette, variable vol. 1-5 ml, Bravo	1.20 / 0.60	0.30 / 0.15	
B10000-1	Capp pipette, variable vol. 1-10 ml, Bravo	3.00 / 0.60	0.60 / 0.15	

### Capp *Bravo* fixed volume (customized volumes available on request)

Cat. No.	Description	Inaccuracy %	Imprecision %
B01-1F	Capp pipette, fixed vol. 1 ul, Bravo	2.50	1.80
B02-1F	Capp pipette, fixed vol. 2 ul, Bravo	2.00	1.00
B05-1F	Capp pipette, fixed vol. 5 ul, Bravo	1.20	0.80
B10-1F	Capp pipette, fixed vol. 10 ul, Bravo	1.00	0.50
B20-1 F	Capp pipette, fixed vol. 20 ul, Bravo	1.00	0.30
B50-1F	Capp pipette, fixed vol. 50 ul, Bravo	1.00	0.30
B100-1F	Capp pipette, fixed vol. 100 ul, Bravo	0.70	0.20
B200-1F	Capp pipette, fixed vol. 200 ul, Bravo	0.60	0.20
B500-1F	Capp pipette, fixed vol. 500 ul, Bravo	0.80	0.15
B1000-1F	Capp pipette, fixed vol. 1000 ul, Bravo	0.80	0.20
B5000-1F	Capp pipette, fixed vol. 5 ml, Bravo	0.60	0.15
B10000-1F	Capp pipette, fixed vol. 10 ml, Bravo	0.60	0.15

### Capp *Bravo* Starter Packs

Cat. No.	Description
SK-03-BB	Three pipettes free of choice, C-05 stand, 3 racks of tips and 1 timer

### Pipette stands

Cat. No.	Description
C-05	Pipette stand for up to 5 pipettes