

Capp *Harmony*

SEROLOGICAL PIPETTES

- DESIGNED FOR ACCURATE AND RAPID FLUID DISPENSING, THESE POLYSTYRENE PIPETTES ARE STERILE, PYROGEN FREE AND COTTON PLUGGED FOR A VARIETY OF LIQUID HANDLING APPLICATIONS.
- THE SHARP, BLACK, HIGHLY LEGIBLE PERMANENT GRADUATIONS PREVENT AMBIGUITY AND GUARANTEE PRECISE READING AND ACCURATE DISPENSING TO + 1-2%. REVERSE GRADUATIONS SHOW VOLUME REMOVED OR REMAINING ON ALL PIPETTES LARGER THAN 1 ML.
- GENEROUS NEGATIVE GRADUATIONS PROVIDE ADDITIONAL PIPETTING VOLUME.
- EACH PIPETTE COMES WITH EFFECTIVE FILTER BARRIER TO PREVENT OVERFLOW OF LIQUID INTO PIPETTE.
- DESIGNED TO FIT ALL AVAILABLE MOTORIZED PIPETTE CONTROLLERS ON THE MARKET.
- PIPETTES ARE SUPPLIED INDIVIDUALLY WRAPPED IN FIBER-FREE COLOUR CODED PAPER AND PLASTIC COMBINATION THAT IS EASY TO OPEN. THE UNIFORM MOUTH PIECE FITS ALL MAJOR PIPETTORS AND IS COLOR CODED FOR EASE OF IDENTIFICATION.

CAPP HARMONY PIPETTES ARE MADE FROM TRANSPARENT MEDICAL GRADE POLYSTYRENE STERILIZED BY GAMMA IRRADIATION. EACH LOT IS TESTED FOR ENDOTOXINS (C.5 EU/ML) (ISO TO993-11). CYTOTOXICITY (ISO 10993-5) AND ARE GUARANTEED NON-HAEMOLYTIC (ISO 10993-4). IDEAL FOR USE WITHIN TISSUE CULTURE APPLICATIONS.



Ordering Info

| CAT.No. | DESCRIPTION | VOLUME | PACKAGE | COLOR |
|---------|---------------------|--------|------------|--------|
| SP-1-C | SEROLOGICAL PIPETTE | 1 ML | 4x200 PCS. | Yellow |
| SP-2-C | SEROLOGICAL PIPETTE | 2 ML | 4x150 PCS. | Green |
| SP-5-C | SEROLOGICAL PIPETTE | 5 ML | 4x100 PCS. | Blue |
| SP-10-C | SEROLOGICAL PIPETTE | 10 ML | 4x100 PCS. | Orange |
| SP-25-C | SEROLOGICAL PIPETTE | 25 ML | 4x 50 PCS. | Red |
| SP-50-C | SEROLOGICAL PIPETTE | 50 ML | 4x 45 PCS. | Purple |