

Compatibility Chart Expell qPCR 96 Well Plate

CAPP Cat. No.	5100219C	5100220C	5100221C	5100319C	5100318C	5100317C
Image						
Standard Well (0.2 mL, height 21 mm)	■	■	■			
Low Profile Well (0.1 mL, height 16 mm)				■	■	■

		qPCR 96 Well Plate						
Bio-Rad		CAPP Cat. No.	5100219C	5100220C	5100221C	5100319C	5100318C	5100317C
Device	Note							
iCycler iQ			■		■			
iQ5			■		■			
MylQ			■		■			
Chromo4			■		■			
Opticon I & II							■	■
CFX96							■	■

		qPCR 96 Well Plate						
Stratagene/Eppendorf/Biometra		CAPP Cat. No.	5100219C	5100220C	5100221C	5100319C	5100318C	5100317C
Device	Note							
Mx4000/3000P/3005P			■		■			
AriaMx							■	■
MasterCycle ep Realplex			■		■			
Optical Thermocycler			■		■			