

EXpell PCR

PCR TUBES & PLATES

CAPP EXpell PCR line comprises of 0.1mL and 0.2mL tubes with plug caps, strips of 8 tubes, strips, plates and optical adhesives. The V-shape bottom allows total sample recovery and the thin-walled wells of polypropylene allow a very good transmission of temperature and minimize condensation.



5100200C



5100210C



5100211C




5100212C



5100213C



5100300C



5100310C



5100311C



5100313C



5100400C

Ordering Information

For standard PCR / qPCR thermocyclers that require regular profile 0.2mL tubes

Cat. No.	Description
----------	-------------

5100200C	ExpellPlus PCR tubes 0.2mL, low-binding, bag, 20x1000 pcs.
5100201C	Expell PCR tubes 0.2mL, bag, 20x1000 pcs.
5100210C	Expell PCR strips of 8 tubes 0.2mL, bag, 20x120 pcs.
5100211C	Expell PCR strips of 8 caps (flat), for 0.1mL and 0.2mL tubes, bag, 30x120pcs.
5100212C	Expell PCR strips of 8 tubes w/ single attached caps 0.2mL, 8x120 pcs.
5100213C	Expell PCR plate of 96 x 0.2mL, 8x 25 pcs., non-skirted
5100214C	Expell PCR plate of 96 x 0.2mL, 8x 25 pcs., frosted, sub-skirted (for ABI cyclers)
5100216C	Expell PCR plate of 96 x 0.1/0.2mL, 8x 25 pcs., clear, fully skirted
5100450C	Expell PCR Shell frame grid 96 x 0.2mL for ABI cycler, 8x24 pcs.
5100451C	Rigid shell frame grid 96 x 0.2mL for ABI cycler, 1 pc

For fast PCR / qPCR thermocyclers that require low profile 0.1mL tubes

Cat. No.	Description
----------	-------------

5100300C	Expell PCR tubes 0.1mL, bag, 20x1000 pcs.
5100310C	Expell PCR strips of 8 tubes 0.1mL, bag, 15x120 pcs.
5100311C	Expell PCR strips of 8 tubes w/ single attached caps 0.1mL, 8x120 pcs.
5100313C	Expell PCR plate of 96 x 0.1mL, 8x 25 pcs., non-skirted
5100314C	Expell PCR plate of 96 x 0.1mL, 8x 25 pcs. white semi-skirted (for Roche LightCycler 480 and COBAS Z 480)
5100315C	Expell PCR plate of 96 x 0.1mL, 8x 25 pcs., frosted, sub-skirted (for ABI fast cyclers)
5100400C	Expell PCR Optical Disposable Adhesive, bag of 100 pcs.