



MICRONIC SAMPLE STORAGE TUBE RANGE

2.50ml Push Cap Tubes

Highest net-working volume on 96-format




Ordering information

Catalog number	Description	Unit
MP32032	2.50ml NC tubes U-bottom in bulk	box/480
MPW32073BC3	2.50ml NC tubes U-bottom in Roborack-96	box/6

* These are Micronic Europe item numbers, for Micronic America item numbers change BC3 to BC

Superior features

- Standard available Non-coded
- Net working volume of 2.50ml at room temperature
- Available in bulk and Roborack-96 (based on the global recognized ANSI / SLAS standards for storage racks)
- 1D rack barcodes are permanently laser etched
- Compatible with the entire range of Micronic TPE and EVA Push Caps
- TPE Push Caps are available in Natural + 12 colors
- Ability to load with standard 8- or 12-channel pipettes (used in combination with extended length pipette tips)
- Class 7 clean room production
- RNase/DNase and Pyrogen free



30 years of experience
in traceable sample
storage

Class 7 clean room
injection molding
and assembly

Production in
the Netherlands and
the United States

To contact Micronic, email to sales@micronic.com or call us at:

Micronic Europe (Europe, Middle-East and Africa): +31 (0)320 277 070 | **Micronic America** (United States, Latin America and Canada): +1 484 480 3372 | **Micronic Asia** (Asia and Pacific): +31 (0)320 277 099