

Disposable PCR Plate 96 well plate – 150 µL



Product description

Riplate® PCR, 96-well, 150 µl, skirted, low profile, clear

Part Number	43001-0305
Well volume	150 µl
Well bottom	V-Bottom
Color well	clear
Color frame	clear
Sterile	no
Re-usable	no
Patient contact	no

Classification

Markets EU

Purity




Produced under cleanroom conditions
(ISO 14644-1 class 8)

Free from DNA, DNase, RNase, and PCR inhibitors.

Materials

Well	PP – Polypropylene
Frame	PC – Polycarbonate

Dimensions	
Plate Length	127.76 mm
Plate Height	15.9 mm
Plate Width	85.47 mm



Packaging	
Box	50 plates
Pallet	84 boxes (4,200 plates)

Specifications are subject to change without prior notice.

In case of any questions, please contact our sales representative via phone, e-mail or our website.

<https://www.ritter-medical.de/en/contact>